

### Product highlights

- Ensures hands-free operation for the staff in case of emergency situations
- Pushing the red TRB button by foot the chair moves into Trendelenburg position



- In Trendelenburg position the entire mattress support platform shall be able to achieve a minimum 12° so that the Patients's head is lower than the circulatory centre point of the body.

#### Technical data

Operating voltage	24V DC
Degree of protection	IPX7
Body material	Thermoplastic
Operating temperature range	- 5C° to + 40C°
Available for	Comfort-3 Eco/Flex/Scale Comfort-4 Eco/Battery/Plus Therapy chairs



Trendelenburg foot-switch  
on both sides of underframe



Trendelenburg foot-switch  
position on Comfort-3  
*installed on both sides of underframe*



Trendelenburg foot-switch  
position on Comfort-4  
*installed on both sides of underframe*

