

Product highlights

- Battery unit ensures chair operation in case of power failure assuring the patient's Comfort and Emergency positions (Trendelenburg, CPR/Shock)
- In case of power failure handset and optional Trendelenburg foot-switches (Art.Nr. 713E005) operates further on from battery
- Battery is automatically charging by control unit. Charging status is indicated by yellow coloured CHARGE LED on the handset.
- Easy to clean plastic surfaces due to compact Battery design
- Proven lead-gel technology
- High IP protection class: IPX6
- Deep-discharge protection (acoustic warning signal when deep discharge level is reached)



Technical data

Capacity	1,2Ah
Output voltage	24V DC
Protection class	IPX6
Operating time	5 hours (tested by one movement in every 10 minutes from seat to TRB and from TRB to seat position, fully charged battery, Patient weight 100kg, room temperature 21°C)
Charging time	approx. 8 hours
Operating temperature	from +10 °C to +40 °C
Available for	Comfort-2 Eco/Flex Comfort-3 Eco/Flex/Scale Comfort-4 Eco Therapy chairs

