

Product highlights

- In case of power failures Battery unit ensures the patient's Trendelenburg & Entry positions
- In case of power failure the Safety Trendelenburg button on control unit and optional Trendelenburg foot-switches (Art.Nr. 713E005) operate further on from battery
- Battery is automatically charging by control unit. Charging status is indicated by yellow coloured CHARGE LED on the control unit.
- Easy to clean plastic surfaces due to compact Battery design
- Proven lead-gel technology
- High IP protection class: IPX6
- Deep-discharge protection (acoustic warning signal when deep discharge level is reached)

Battery unit
lead-gel technology



Black Entry button
to set sitting position in
case of power failure



Safety TRB button
operates further from
Battery



Yellow CHARGE LED
indicates battery charging
on Control unit

Technical data

Capacity	1,2Ah
Output voltage	24V DC
Protection class	IPX6
Operating time	5 hours (tested by one movement in every 10 minutes from seat to TRB and from TRB to seat position, fully charged battery, Patient weight 100kg, room temperature 21°C)
Charging time	approx. 8 hours
Operating temperature	from +10 °C to +40 °C
Available for	Comfort-4 Plus Therapy chairs

