

Product highlights

- Battery unit ensures chair operation in case of power failure assuring the patient's Comfort and Emergency positions (Trendelenburg, CPR/Shock)
- In case of power failure handset and optional Trendelenburg foot-switches (Art.Nr. 713E005) operates further on from battery
- Battery is automatically charging by control unit. Charging status is indicated by yellow coloured CHARGE LED on the handset.
- Easy to clean plastic surfaces due to compact Battery design
- Proven lead-gel technology
- High IP protection class: IPX6
- Deep-discharge protection (acoustic warning signal when deep discharge level is reached)



Technical data

Capacity	2,4Ah
Output voltage	24V DC
Protection class	IPX6
Operating time	10 hours (tested by one movement in every 10 minutes from seat to TRB and from TRB to seat position, fully charged battery, Patient weight 100kg, room temperature 21°C)
Charging time	approx. 16 hours
Operating temperature	from +10 °C to +40 °C
Available for	Comfort-3 Eco Comfort-4 Eco/Battery Therapy chairs

