

### Product highlights

- “Z” shape foldable side rail is designed for maximum Patient Safety – effective protection for patient from falling off the chair
- By pressing the red handle easily fold the side rail down for patient getting in and out from the chair
- Fixed on both sides of seat frame. Optionally it can be also mounted on DIN rail (with DIN rail Art.Nr. 713M063) installed on both sides of seat.
- No sharp edges as additional safety feature
- High durability plastic injection arms
- High Strength structure



Ergonomic handle design  
for Easy Lock & Unlock

#### Technical data

Load capacity	15 kg
Size of rail (length x height in fold up/ fold down position)	120 x 38 cm / 120 x 11 cm
Material	Polyamid 6 + Glassfiber %30 + anodized aluminium + powder coated steel profile
Available for	Comfort-2 Eco/Flex Comfort-3 Eco/Flex/Scale Comfort-4 Eco/Battery/Plus Therapy chairs

Foldable side rail  
on both sides of seat



**UP position**

*Folded up position ensures effective protection for patient from falling off the chair during treatment or transport*



**DOWN position**

*Lowering of side rail allows the patient to easily enter or leave the chair and free access to the patient by staff*



