

Product highlights

- “Z” shape foldable side rail is designed for maximum Patient Safety – effective protection for patient from falling off the chair
- By pressing the red handle easily fold the side rail down for patient getting in and out from the chair
- Fixed on both sides of seat frame. Optionally it can be also mounted on DIN rail (with DIN rail Art.Nr. 713M063) installed on both sides of seat.
- No sharp edges as additional safety feature
- High durability plastic injection arms
- High Strength structure



Ergonomic handle design
for Easy Lock & Unlock

Technical data

Load capacity	30 kg
Size of rail (length x height x width)	92,5 x 34,5 cm x 3,6 cm
Material	Reinforced Polypropylene + powder coated steel frame
Available for	Comfort-2 Eco/Flex Comfort-3 Eco/Flex/Scale Comfort-4 Eco/Battery Therapy chairs

Tiltable side rail
on both sides of seat



UP position

Tilted up position ensures effective protection for patient from falling off the chair during treatment or transport



DOWN position

Lowering of side rail allows the patient to easily enter or leave the chair and free access to the patient by staff



