

Product highlights

- Plastic side rail set is designed for maximum Patient Safety – effective protection for patient from falling off the chair
- By pulling the red handle easily tilt the side rail down for patient getting in and out from the chair either from front or side
- High durability blow molded, reinforced plastic
- Flame retardant, antibacterial material
- No sharp edges as additional safety feature
- Crush resistant / High Strength structure
- Fixed on both sides of seat frame
- Smooth surface – easy to clean



Ergonomic handle design
for Easy Lock & Unlock

Tiltable side rail
on both sides of seat



Technical data

Load capacity	30 kg
Size of rail (length x height x width)	92,5 x 34,5 cm x 3,6 cm
Material	Reinforced Polypropylene + powder coated steel frame
Available for	Comfort-2 Eco/Flex Comfort-3 Eco/Flex/Scale Comfort-4 Eco/Battery Therapy chairs

UP position

Tilted up position ensures effective protection for patient from falling off the chair during treatment or transport



DOWN position

Lowering of side rail allows the patient to easily enter or leave the chair and free access to the patient by staff



